

PRODUCT NEWS


22

THE ARNOUDIZER



33-54730

Arnoudizer, incl. fixation screw
Arnoudizer, inkl. Fixationsschraube



An excellent device for molar distalization. The upper molars will be derotated, uprighted and distalized in only two or three months. By converting the system in a transpalatal bar, it can be used as an anchor for the molars. The anterior teeth can be retracted and aligned simply.

Eine hervorragende Vorrichtung zur Molarendistalisierung. Die oberen Molaren werden in nur zwei oder drei Monaten derotiert, aufgerichtet und distalisiert. Durch die Umwandlung des Systems in einen transpalatinalen Bogen, kann der Arnoudizer als Verankerung für die Molaren verwendet werden. Die Prämolaren können einfach distalisiert und ausgerichtet werden.

BENEFIT® .KFO

Material: Stainless steel
Packaging: 1ea

Material: Rostfreier Stahl
VPE: 1/ Pkg

Developed in cooperation with:
Dr. Arnoud Schreurs, Orthodontist,
Genk, Belgium

Entwickelt in Zusammenarbeit mit:
Dr. Arnoud Schreurs, Orthodontist,
Genk, Belgium

Manufacturer | Hersteller:

psm MEDICAL SOLUTIONS
Moltkestraße 41
78532 Tuttlingen, Germany
Phone +49 (74 61) 9 66 37-0
Fax +49 (74 61) 9 66 37-29
E-Mail info@psm.ms
Web www.psm.ms